## ergotron® moving you forward

## Neo-Flex<sup>®</sup> Dual Monitor & Keyboard Wall Mount

- Combine with the Neo-Flex Keyboard Wall Mount (45-403-062) to create a complete imaging workstation
- Easily access PACS data
- No more light boxes—work directly with digital records
- Low profile makes it perfect for space-constrained areas
- Dual-monitor configuration increases
   productivity
- Comfortable viewing: Tilt screens to reduce glare and improve sight lines
- Includes 5-year warranty



## MAXIMIZE SPACE IN IMAGING ENVIRONMENTS!

This low-profile digital imaging solution enhances PACS implementation by providing convenient access to patient data.

The slide-out, low-profile keyboard solution easily folds up out of the way when not in use. The back-tilt tray enhances user comfort, while its side-to-side adjustment capability provides positioning flexibility to accommodate work space requirements.

## Neo-Flex<sup>®</sup> Dual Monitor & Keyboard Wall Mount

- Create a complete digital imaging station by combining the Neo-Flex Keyboard Wall Mount with the Neo-Flex Dual Monitor Wall Mount.
- For a single monitor solution, pair the Keyboard Mount with the Neo-Flex Wall Mount Lift (60-577-195).



FOLDS TO 5"





TRAY CAN EXTEND 5" (13 CM), PLACING THE KEYBOARD UP TO 19.25" (49 CM) FROM THE WALL



Part # (color)	2
Typical LCD Size	

45-403-062 (white) ≤ 24″

LCD SIZE APPROXIMATION. DIMENSION IS DIAGONAL MEASUREMENT OF SCREEN. LCD SIZE MAY BE EXCEEDED AS LONG AS SCREEN WEIGHT DOES NOT EXCEED MAXIMUM WEIGHT CAPACITY OF THE MOUNT.

CALACITY OF THE MOONT.	
Weight Capacity	≤ 50 lbs (22,7 kg) 25 lbs (11,3 kg) per pivot Keyboard ≤ 5 lbs (2,3 kg) See website for details
Lift	—
Tilt	LCDs 180° Keyboard 90°
Pan	Keyboard 120°
Rotation	LCDS 360° P/L
VESA	MIS-D
Depth	Keyboard 5" (13 cm)
Extension	Keyboard 19.25" (49 cm)
Shipping Dimensions	28-514-800: 30" x 6" x 6" (76,2 x 15,2 x 15,2 cm) 45-403-062: 20.75" x 15.5" x 5.7" (53 x 39,5 x 14 cm)
Shipping Weight	28-514-800: 12 lbs (5,5 kg) 45-403-062: 13 lbs (5,9 kg)
Warranty	5 years

 Download additional resources at ergotron.com.
 © Ergotron, Inc. 07.18.2019

 For more information:
 Literature made in the U.S.A

 NORTH AMERICA: 800.888.8458 / +1.651.681.7600 / sales@ergotron.com
 Content subject to change.

 EMEA: +31.33.45.45.600 / info.eu@ergotron.com
 Ergotron devices are not in

 APAC: apaccustomerservice@ergotron.com
 Patent information available

 OEM: info@oem.ergotron.com
 Etgotron.com

© Ergotron, Inc. 07.18.2019 Literature made in the U.S.A. Content subject to change. Ergotron devices are not intended to cure, treat, mitigate or prevent any disease. Patent information available at www.ergotron.com/patents.

